




The DRY FLOATATION® BARIATRIC KIT

by ROHO, Inc.

The DRY FLOATATION Bariatric Kit, which works in conjunction with the ROHO DRY FLOATATION Mattress System, has been designed to address the needs of the clinician/caregiver and the compromised obese patient by creating an environment in which the patient can comfortably lie or sit for extended periods of time without the fear of compromised blood flow and the resulting effect of tissue damage/breakdown. This ROHO mattress system achieves this by allowing your bariatric patients to benefit from compression therapy which addresses edema associated with chronic wounds and flaps, and by a system which conforms to your patient's body without deforming the vascular bed.

The DRY FLOATATION Mattress System for your bariatric patients - the **ORIGINAL non-powered, adjustable, zoned, pressure-reducing air mattress overlay**. The Mattress System where weight doesn't matter.

1-800-851-3449 ♦ 618-277-9173 ♦ Fax 618-277-9561 ♦ www.crownthera.com



ROHO MATTRESS
SYSTEM



HIGH PROFILE®
QUADTRO®

CROWN®
THERAPEUTICS, INC.
The Promise of Discovery™



ENHANCER®




HIGH PROFILE®
CUSHION

Features and Benefits

Zoned: The base of the ROHO Mattress System comes in four sections. Add on sections of 3 x 10 and 4 x 10 cells allow the bariatric patient the therapeutic benefits of ROHO DRY FLOATATION technology. Each section can then be adjusted to address the specific needs of your patient's body shape and size.

Adjustable: The three (3) inch high air cells can be adjusted to shape and conform to your patient's body and provide a total contact environment ideal for facilitating wound healing.

Ventilation: Individual air cells allow for air circulation between the cells.



Drainage Holes: To pull fluids away from the highly incontinent patient, and to address the problems of maceration associated with wound healing.

Durable Neoprene Rubber: Flame resistant and resistant to urine, oils and medicines. Easily cleaned with soap or water and can be steam autoclaved. Product can be repaired extending the overall life of the cushion.

Lightweight: The ROHO Bariatric Kit weighs approximately one (1) pound per section or five (5) pounds for a complete system*.

*Weight based on a 3 x 10 cell section.

24-Month Limited Warranty

Specifications

The ROHO DRY FLOATATION mattress comes complete with hand inflation pump, repair kit and operating instructions.

Construction—Flame resistant neoprene rubber

Weight—Approximately one (1) pound per section

Size—Systems are available to fit 39" x 80", 48" x 80" and 54" x 80" beds.

Cleaning—Clean with soap and water, or can be steam autoclaved

There is no weight limit for this mattress system.

Applications

Can be used on patients with any stage ischemic ulcer.

CROWN[®]
THERAPEUTICS, INC.

100 Florida Avenue • Belleville, IL 62221-5429 USA